


## Mouth Orthodontic Expanders Adults Plate For Dental Alignment

Our Product Introduction

for more products please visit us on [adschinadentallab.com](http://adschinadentallab.com)

### Basic Information



### Product Specification

- Turnaround Time: 3 Days In Lab
- Warranty: 5 Years
- Shipping Way: UPS (preferred), Fedex And DHL Linex. Armex Are Also Available
- Payment: T/T, HSBC, PayPal, MoneyGram, Western Union
- Highlight: Dental orthodontic expander adults, Mouth orthodontic expander adults, Dental orthodontic expander plate



### More Images



## Product Description

### Description

Orthodontic Expanders: An Advanced Solution for Dental Alignment

As a manufacturer located in China, we are proud to offer our latest innovation to merchants worldwide – our Orthodontic Expanders. Our expanders are designed to provide an advanced solution for dental alignment, offering a wide range of benefits to both patients and practitioners alike.

Our Orthodontic Expanders are manufactured using high-quality materials and advanced technologies, ensuring they are both effective and durable. We take great care in producing our expanders to meet the highest standards of quality, ensuring they are safe to use and can be relied upon for optimal performance.

One of the key advantages of our Orthodontic Expanders is their ability to treat a range of orthodontic issues. The expanders can be used to treat an array of dental conditions, such as crowded teeth, crossbite, and overbite, among others. This versatility makes the expanders an attractive option for dentists seeking to address a range of patient needs.

Our Orthodontic Expanders are also designed with patient comfort in mind. We know that dental procedures can be stressful and anxiety-provoking for patients, and we have taken steps to ensure that our expanders are as comfortable as possible. With minimal discomfort or pain, patients can enjoy improved dental alignment and other associated benefits.

Another benefit of our Orthodontic Expanders is their non-invasive nature. With no need for invasive surgical procedures, patients can enjoy a swift and smooth recovery with minimal discomfort or downtime. This non-invasive nature also makes the process more affordable and less time-consuming, allowing patients to undergo treatment without undue burden.

Overall, our Orthodontic Expanders offer a range of benefits to patients and practitioners alike, making them an ideal solution for the treatment of orthodontic conditions. With their advanced design, high-quality materials, and proven effectiveness, our expanders are sure to provide an excellent return on investment for merchants seeking to offer advanced dental solutions to their customers.

If you are interested in learning more about our Orthodontic Expanders or ordering a batch for your business, please do not hesitate to contact us. We are always happy to answer any questions you may have and offer advice and support wherever needed. Thank you for considering our product – we look forward to working with you to improve dental health and alignment worldwide.

### Production Introduction:

A palate expander is a custom-made dental device. It's used to widen the upper jaw, creating more space for teeth. Dentists or orthodontists may use a palate expander before placing dental braces.

The expander has two halves that attach to the upper teeth on each side. The two halves are joined together by a "suture" in the center of the roof of the mouth. The RPE is attached to the upper back teeth and eases the suture apart, which makes the jaw wider. As the jaw expands, new bone fills in between the two halves of the palate.

### Product Parameter

Customized

### Product Feature and Application

An orthodontist may use a palate expander to widen the roof of the mouth if someone has a narrow upper jaw. This creates more room for the teeth.<sup>2</sup>

Dentists recommend palate expanders for people with:

- Malocclusion
- Dental crowding

### Types of Palate Expanders

There are several types of palatal expanders, including:

Removable palatal expander  
Hyrax rapid palatal expander  
Haas expander  
Quad-helix appliance

### Deliver, Shipping and Serving

We accept STL files.

We accept normal impressions, you can send us by Fedex, UPS, DHL or Linex

We finish the jobs in 4-5 days once we receive your cases

We accept TT, Paypal

**Why choose us?**

We offer a wide selection of Orthodontic Appliances to meet all your needs.

The purpose of our company is to create wealth for society, profit for customers, value for the company, and opportunities for employees.

We are committed to providing transparent pricing and clear communication about our products and services.

We are a Orthodontic Expanders manufacturer integrating R&D, manufacturing, sales and after-sales service. In the production process, we strictly implement production standards to ensure product quality.

Our products are backed by a comprehensive warranty to ensure your complete satisfaction.  
In recent years, our company has been strengthening the technical innovation of equipment and process transformation for the purpose of energy saving, consumption reduction and emission reduction.  
We take pride in our ability to provide fast and efficient service, without sacrificing quality or attention to detail.  
We actively explore and innovate, adhere to strict delivery time limit, reasonable price standard, excellent product quality and perfect service.  
Our products are of professional quality and come with excellent customer service.  
Over the years, we have been implementing the strategic goal of becoming a stronger and better enterprise with our brand advantages and advanced management experience.  
Hot Tags: orthodontic expanders, China orthodontic expanders manufacturers, suppliers, factory, palate expander lower jaw, bonded palatal expander, Frankel Dental Devices, palate expander and braces at same time, attachment partial denture, mouth expander for braces



+8613430498390

info@chinaadsdentallab.com

adschinadentallab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103